User Manual
Fingertip Pulse Oximeter
Mode d'emploi
Oxymètre de pouls (doigt)
Benutzerhandbuch
Fingerspitzen-Pulsoximeter
Gebruikshandleiding
Vingerpit Pulseoximeter

Version:Ver.1.0

Printed Version Number: ENSM80-84201909161

Applicable Models: YK-80A, YK-80B, YK-80C,

YK-81A, YK-81B, YK-81C, YK-82A, YK-82B, YK-82C,


YK-83A, YK-83B, YK-83C, YK-84A, YK-84B, YK-84C

Xuzhou Yongkang Electronic Science Technology Co., Ltd.
4F Building C8, 40 Jingshan Road, Economic and Technological Development Zone, 221000 Xuzhou, PEOPLE'S REPUBLIC OF CHINA

EC REP Company Name: Proline GmbH
Company Address: Brühnstr.56,40239,Duesseldorf Germany
Competent Authority: Code DE/CA200
Fax : 0049 2131 4051968-9
Tel:0049 2131 4051968-0
Email : med@eulinis.eu

2. In menu 2, When the "*" signal is shown on the "Sounds Setup", press the button (>3s) and enter into the menu 2. Press the button ('<1s) can select item,then press the button (>3s) to setup data,choose "+or"- to plus or minus values.

Note: If the matching refers to clamp oximeter is basic, not equipped with parameter setting function

Detailed descriptions of product functions:

- Display Type: OLED display
- SpO2: Measurement range: 70% ~ 99% Accuracy: 80% ~ 99% ±2% (Including 80%) ; 70% ~ 79% ±3% ; Below 70% no requirement;
- Resolution: 1PPM
- PR: Measurement range: 30BPM ~ 240BPM Accuracy: ±1BPM or ±1% (the larger one)
- Parameters of LED probe

wavelength	radiation power
RED	660±2nm 1.8mW
IR	905±10nm 2.0mW

- This parameter can be especially useful to clinicians.
- Pt index: measurement scope: 0~20 (optional)
- Power: two AAA 1.5V alkaline batteries
- Automatic standby: the product shuts off by itself when no finger is in the product about 8 seconds
- Automatic startup: every 5 s instrument will automatically detect the signal, after the hole with my finger, timely automatically boot;(optional)
- Accelerometer function: finger movement, the screen display will change with the accelerometer changes (optional)
- Dimension: 58mm×36mm×35mm

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Ver.1.0

Fingertip Pulse Oximeter

This Fingertip Pulse Oximeter is a kind of innovated medical device with non-invasive and continuous features for artery SpO2 and PR detection. Being portable, it is able to measure SpO2 and PR values quickly and precisely.


General Description

Haemoglobin Saturation is the percentage between the capacity of Oxyhemoglobin (HbO2) that compounded with oxygen and that of all combinatory haemoglobin (Hb) (obin HbO2) in blood. In other words, it is the saturation of Oxyhemoglobin in blood. It is a very important physiological parameter for Respiratory and Circulatory system. It is also a key indicator of oxygen supply and oxygen saturation in human blood. Moreover, factors such as Automatic Organic Regulation Malfunction caused by anaesthesia, trauma resulted from major operation and some medical examination can also cause problems in oxygen supply, which might reduce human haemoglobin saturation. As a result, such symptoms as megrim, vomiting and asthenia might appear to patients. Hence, it is very important to know hemoglobin saturation of patient timely in clinical medical aspects. The fingertip pulse oximeter features in small volume, low power consumption, convenient operation and portability. It is only necessary for patient to put one finger into fingertip photoelectric sensor for measurement, which will directly show measured value of hemoglobin saturation. It has been proved in clinical experiments that it possesses rather high precision and repeatability.

Measurement principle

The principle of the oximeter is as follows: An experience formula of data process is established by exerting Lambert Beer Law according to Spectrum Absorption Characteristics of Reductive hemoglobin (R.Hb) and oxyhemoglobin (O2 Hb) in glow and near-infrared zones. Operation principle of the instrument is to combine Photoelectric Oxymeter Inspection Technology with Capacity Pulse Scanning and Recording Technology, so that two lights with different wavelength (660nm glow and 940nm near infrared light) can be focused onto the same sensor site simultaneously. Then the measured signal will be obtained by a photosensitive element, information acquired through which will be shown on two groups of LEDs through process in electronic circuits and microprocessor.

Diagram of Operation Principle

Precautions for use

- Do not use the Fingertip Pulse Oximeter together with MRI or CT equipment.
- Excessive hazard: Do not use the Fingertip Pulse Oximeter in an explosive atmosphere.
- The Fingertip Pulse Oximeter is intended only as an adjunct in patient assessment. Doctors should make diagnosis in conjunction with clinical manifestation and symptoms.
- Check the Fingertip Pulse Oximeter sensor application site frequently to make sure that the circulation and skin integrity of patient are under good condition.
- Do not stretch the adhesive tape while applying the Fingertip Pulse Oximeter sensor. This may cause inaccurate reading or skin blisters.
- Please read the manual carefully before your operation.

- Do not stretch the adhesive tape while applying the Fingertip Pulse Oximeter sensor. This may cause inaccurate reading or skin blisters.
- Please read the manual carefully before your operation.
- Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status, and correct alignment at least every 2 hours.

- Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid.

- Significant levels of dysfunctional hemoglobin (such as carboxyhemoglobin or methemoglobin) may cause inaccurate reading.

- Intravascular dyes such as indocyanine green or methylene blue may cause inaccurate reading.
- Some sensors may be adversely affected in the presence of high ambient light. Please shield the sensor area (with a surgical towel or direct sunlight, for example) if it is necessary.

- Unexpected action may cause inaccurate reading.

- Medical signal with high frequency or interference caused by defibrillator may lead to inaccurate reading.

- Venous pulsations may cause inaccurate reading.

Product Accessories

- One hand lace
- One user manual

Guidance and manufacturer's declaration-electromagnetic radiation-for other EQUIPMENTS and SYSTEMS

The Fingertip Pulse Oximeter is designed to be used in specified electromagnetic environment. Users of the Fingertip Pulse Oximeter must use it in the following environments.

Classification

1. Management Class for Medical Devices:

II equipment

2. Anti-electric Shock Type: Internally powered equipment

3. Anti-electric Shock Degree: Type BF equipment

Maintenance and Preservation

1. Replace the batteries timely when low voltage lamp is on.

2. Clean the surface of fingertip pulse oximeter before it is used to measure patients.

3. Replace the batteries inside if you will not operate the Oximeter for a long time.

4. It would be better to preserve the product in -10~40°C (-14~104°F) and humidity is 10%-80%.

5. It is recommended that the product should be kept dry anytime. A wet ambient might affect its lifetime and even damage the product.

6. Please follow the law of the local government to deal with used batteries.

Symbols and Definitions

	BF type application part	IP22	IP degree	SN	Serial number
	Separate collection		Humidity range		Date of manufacture
	Reference manual		Keep dry		Manufacturer
	Cautions		Product certification		European union representative
	Standby		Avoid sunlight		
Up toward					

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- 1. Mode de mesure : 70% ~ 99% Accuracy : ±2% (including 80%) ; 70% ~ 79% ±3% ; Below 70% no requirement;
- 2. Mode de fréquence : 30BPM ~ 240BPM Accuracy : ±1BPM or ±1% (the larger one)
- 3. Paramètres de LED probe

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