

**VINNO E35**  
Affordable Luxury



**VINNO**

**VINNO Technology (Suzhou) Co., Ltd.**

5F, A Building, No.27 Xinfu Rd, Suzhou Industrial Park, 215123, China

Tel: +86 512 62873806

Fax: +86 512 62873801

Email: [vinno@vinno.com](mailto:vinno@vinno.com)

Website: [www.vinno.com](http://www.vinno.com)

VINNO reserves the rights to revise the technical specifications if needed.

**VINNO**

**VINNO E35 is an affordable luxury ultrasound system for general purpose with delicate design. Its versatility meets the demands of any clinical setting.**



- Superb Imaging
- Delicate Design
- Agile Workflow
- Expanded Application



# Superb Imaging

## Innovative RF Platform

Due to VINNO's innovative RF Platform, the system can acquire 40 times amount of raw signal for back-end process which represents better resolution and powerful post process capability.



Precise IMT Measurement within 0.1mm



Low Velocity Color Flow by RF Platform

## Advanced Image Processing Technologies

**VFusion** is spatial compound technology from VINNO which helps to improve the image information, enhance contrast resolution and reduce the noises.

**VSpeckle** is an intelligent algorithm which can remove noise, enhance edge and contrast of the image.



# Delicate Design




Make your working environment more sophisticated and compact design.



## Expanded Application

### Easy Compare




Provide compared information of last exam to make confident diagnosis

### 4D-HQ



HQ Light rendering generates realistic images of the fetal structure for increased diagnostic confidence

### Auto Follicle




Just one click, the numbers and sizes of follicles can be calculated automatically




## Agile Workflow


Smart touch panel operation and self-defined exam protocol makes your daily work high efficient.




Self-defined preset



Self-defined measurement items




Fast review by sliding on touch screen




Built-in report template


# Image Gallery




Pregnancy Yolk Sac Formation




Fetal Heart




Carotid Plaque




Carotid




Hepatic Vein




Blood Flow Around Thyroid Nodule



Pericardial Effusion




Aortic Regurgitation Flow




Follicle Of Ovary



Fetal Mouth



Adenoma



Fetal Face